CIP-Anlage / CIP Plant

CIP-Anlagen

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Vorderansicht LetzAnalyze-CIP /
Front View LetzAnalyze-CIP

CIP installations

The CIP systems are individually designed and conceived with regard to differing cleaning tasks. Usually, the CIP installations comprise one reservoir vessel, various dosing stations, a hygienic pump, a heat exchanger to heat the media, and the required measurement sensors, such as flow, conductivity, pH-value etc. Depending on the application, the cleaning medium is discarded or completely / partially re-circulated again after cleaning in order to reduce the amount discarded. Identical to the water treatment installations, specific design features such as the use of 316L stainless steel, the reduction of dead spaces, the installation of aseptic connections and the maintenance of a gradient are implemented in order to prevent the growth of microbes and to simplify cleaning of the piping that comes in contact with the product.

Features: LETZNER CIP-Installations

- Complete residue emptying
- Program freely configurable
- Hygienic design

LetzAnalyze-CIP

A validation is essential in order to achieve reproducible, successful cleaning with a CIP cleaning procedure. The analysis of the final rinse is an opportunity for sampling to be done in the context of the cleaning validation for the CIP cleaning procedure and, later, in routine monitoring.

LetzAnalyze-CIP is a measurement station that can be used to simultaneously check different CIP systems, in that the final rinse is analyzed fully automatically. The measurement station consists of the following, main components:

- Online conductivity probes (one probe per CIP system)
- Online TOC measurement device LetzTOC
- Fully automated sampling for offline analysis in the laboratory
- Fill level- and mounting monitoring of the bottles taken for sampling
- Installation of a hygienic pump for non-circulating CIP systems
- Connections, piping, and programming for fully automated Pure Steam sterilization of the system
- Connections, piping and programming for fully automated drying of the system with purified compressed air after Pure Steam sterilization
- PLC, data documentation and remote maintenance with industrial PC (21 CFR Part 11)

The measurement station can be used as a mobile unit. The installation serves to accelerate the cleaning validation and to continually monitor the complete CIP process, corresponding to the PAT (Process Analytical Technologies) initiative of the FDA. The cleaning validation for new CIP processes is accelerated, because the online TOC analysis eliminates lab analysis times for TOC. A reduction in the use of cleaning agents and a greater potential for saving water are possible. The fully automated sampling that is integrated into the measurement station can clearly reduce the risk of contamination during sampling. The complete sampling unit can be mounted and dismantled via aseptic clamp connections. The entire unit can thereby be brought into the lab for analysis without removing the bottles individually. Furthermore, the system can be sterilized with Pure Steam before the next cleaning procedure in order to inhibit microbial growth in the measurement station.