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## Standby operation

If there is no demand for WFI (no tank request), the temperature in the pressure columns is lowered to  $> 90\text{ }^{\circ}\text{C}$ . By doing this, the microbiological quality is assured even in standby phases and new WFI can be produced at the required temperature as quickly as possible.

## Production of Water for Injection by membrane processes

A reverse osmosis membrane stage with downstream electro-deionization and ultrafiltration form the crucial preparation steps. The ultrafiltration is used here as a security barrier to retain microorganisms and endotoxins. Since this is, unlike in the preparation by distillation, to be a cold manufacturing process, ensuring compliance with the microbial quality requirements plays a major role. For this reason influences which support the growth of microorganisms should be reduced and measures which prevent the growth should be integrated. Letzner meets the high quality requirements with a comprehensive assessment of the plant design and additional qualification measures.

Vollisolierte Destillationsanlage /  
Completely Insulated Distillation Unit

### Features: LETZNER WFI production by membrane processes

- The installation of a UV lamp with a wavelength of 254 nm ensures a significant reduction of microorganisms present in the feed water. A second UV lamp with a wavelength of 184 nm prevents indirectly the microbial growth by oxidizing organic substances to  $\text{CO}_2$ . Thus the organic substances are no longer available as a nutrient source
- In standby phases, starting with the pre-treatment, water is constantly circulating through all pipes in contact with the product, so that stagnant water in which microorganisms can multiply is excluded
- The number of connections is reduced to a minimum. At locations where connections are required, aseptic connections are used exclusively
- Only diaphragm valves are used
- For the ultrafiltration ceramic modules are used
- By using 316L stainless steel, the entire system can be sanitized regularly with hot water  $> 80\text{ }^{\circ}\text{C}$ . In addition, certain parts of the system can be charged with ozone (e.g. ultrafiltration)