## Online TOC-Messung

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
### Online TOC Measurement

The analytical procedure with the LetzTOC measurement device allows the entire organic carbon content in a sample to be determined. The system was designed and developed especially for monitoring pharma-water (PW, WFI and Pure Steam condensate). All requirements of the European, United States, and Japanese pharmacopoeias have been fulfilled. The analyzer can measure water from two different sampling locations online. In the first step the inorganic carbon (TIC) is removed from the water sample and is measured. Using UV-light, the organic substances contained in the water are then oxidized to CO<sub>2</sub>, which is removed from the water and is led to the detector. The detection takes place through irradiation with infrared light, which is absorbed by the CO<sub>2</sub> molecules. The diminishment of light is recorded by the detector and is directly proportional to the TOC concentration. Based on the wide measurement range, the analysis of somewhat more contaminated samples is also possible. Thus, for example, LetzTOC can also be used for online monitoring of the final rinse with CIP-systems by which higher TOC values can arise, for example as a result of contamination.

The verification of organic substances by means of infrared light reduces the risk of determining false TOC values. This is because conductivity sensors, which are primarily used for the online TOC measurement, record all ions contained in water and not only the CO<sub>2</sub> resulting from UV-oxidation, respectively the hydrogen carbonate (HCO<sub>3</sub><sup>-</sup>) and carbonate (CO<sub>3</sub><sup>2--</sup>) ions formed by dissociation. The infrared wavelength used for the detection in the LetzTOC measurement device is selective for CO<sub>2</sub>. In addition, the CO<sub>2</sub> is removed from the water before the analysis, whereby, in the event of contamination, other substances do not even arrive at the detector.

### Features: LetzTOC

- TOC determination by infrared spectroscopy - NDIR detection
- Sequential measurement of two independent sampling points
- Measurement of cold and hot water
- Display and trend of the TOC value
- Simple data analysis via USB and archiving with the evaluation program LetzData
- Certified test solutions for Calibration and System Suitability Test (SST)

Visualisierung TOC-Messgerät / Visualization of the TOC Measurement Device